

<b>Product</b>	Name	Organoid® BAMBOART MARGERITTA medium on self-adhesive backing
	Article number	BMBMAR0045SKTRW
	<b>Producer</b>	Address ORGANOID TECHNOLOGIES GMBH Nesselgarten 422 / Top 5 A-6500 Fließ Österreich
<b>Description and applications</b>	Name	Organoid® BAMBOART MARGERITTA medium on self-adhesive backing
	Description	Natural plant material in translucent density on translucent self-adhesive foil, produced climate-neutral in Tyrol, Austria
	Applications	Suitable for applications where the product is visible from both sides. Direct application on smooth surfaces such as glass, metal or composite materials
<b>Plant material and density</b>	Plant material	Grey lichen and approx. 8 marguerites/m <sup>2</sup>
	Density	Medium, not full coverage, self-adhesive backing is part of the design
<b>Size and weight</b>	Size	Rolled goods, width 1 360 mm
	Thickness	0.8 mm
	Weight	400 g/m <sup>2</sup>
<b>Product requirements</b>	Food-safety	Binding agent safe for usage in contact with food in accordance with EU-regulation 1935/2004
<b>Backing material</b>	Material Properties	Acrylic adhesive PET Translucent, strong adhesive Poses no risks for human health and environment PVC-free Contains no dangerous substances in accordance with annex 2, § 3 of the regulation (EC) 1907/2006 heat-resistant up to 120° C
<b>Technical implementation</b>	Cut	Use scissors, a cutter or a plotter
	Edges	Allow for imperfections 1-2 mm around the edges
	How to install	Draw the release paper from the self-adhesive foil, place on chosen surface (dry application)
	Pay attention to	Exact placement of the foil is important, as readjustments are not possible due to the strong adhesive qualities of the product
		Acclimatisation of Organoid natural surfaces: - Bring product to room temperature before

unpacking

- Store product protected from light, at about 20°C and 50-65 % relative humidity

Processing guidelines  
available at:

<https://organoid.corestad.com/verarbeitung/>  
(= *current link*)